

57. 15-05-2009 08:43:45

Roatan


Information from dive computer

|                |                            |                  |                           |                   |               |
|----------------|----------------------------|------------------|---------------------------|-------------------|---------------|
| Dive no.       | <b>57</b>                  | DC Dive no.      | <b>31</b>                 | Sampling Rate     | <b>1 sec</b>  |
| Start time     | <b>15-05-2009 08:43:45</b> | Surface Interval | <b>01:59:00</b>           | Altitude mode     | <b>A0</b>     |
| End time       | <b>15-05-2009 09:37:45</b> | ID no.           | <b>72101106</b>           | Personal settings | <b>P0</b>     |
| Dive time      | <b>00:54:00</b>            | Model            | <b>D6</b>                 | Dive type         | <b>Nitrox</b> |
| Max depth      | <b>13.5 m</b>              | Personal info    | <b>knuspar@knuspar.dk</b> | Alarms Recorded   | <b>1</b>      |
| Oxygen Percent | <b>33% 40%</b>             |                  |                           | OLF               | <b>12 %</b>   |


|             |   |               |              |             |                        |
|-------------|---|---------------|--------------|-------------|------------------------|
| Weather     | <b>Clear Sky / Night</b>                    | Temperatures  |              | Buddy       | <b>Thomas Handberg</b> |
| Visibility  | <b>~20 M</b>                                | Start of dive | <b>29 °C</b> | Boat name   | <b>Cocoview IV</b>     |
| Weight used | <b>7 kg</b>                                 | Max depth     | <b>28 °C</b> | Dive Master | <b>Kirk</b>            |
| Dive gear   |   | End of dive   | <b>28 °C</b> | Custom 1    |                        |
|             | 5mm Shorty + rental, Canon G10 + UW Housing |               |              | Custom 2    |                        |
|             |   |               |              | Custom 3    |                        |
|             |   |               |              | Custom 4    |                        |
|             |   |               |              | Custom 5    |                        |

Cylinders


| Cylinder size    | Oxygen Percent | Start        | End          | Average depth | SAC Rate     | Gas used Pressure used | Gas used Rate used | Gas used Volume gas used |
|------------------|----------------|--------------|--------------|---------------|--------------|------------------------|--------------------|--------------------------|
| <b>12 liters</b> | <b>33%</b>     | <b>0 bar</b> | <b>0 bar</b> | <b>9,1 m</b>  | <b>0 SLM</b> | <b>0 bar</b>           | <b>0 bar/min</b>   | <b>0 liters</b>          |



Tissue saturation on start



Tissue saturation on end



Notes